Tutoriel Study-card Shanoir

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1. Create Center
2. Create Investigator
3. Create Study
4. Add the new Center-Investigator to the Study
5. Create Acquisition Equipement (for the center)
 - a. Create MrModel
 - b. Create Manufacturer
 - c. Create Acquisition Equipement
6. Create Coil (for the acquisition equipement)

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Study card field	Possible values
Protocol name:	...
Protocol Comment:	...
Dataset modality type:	MR dataset
Processed dataset type:	Reconstructed Dataset
Cardinality of related subjects:	Single-Subject Dataset
MR Sequence k-space fill:	Conventionel cartésien
Transmitting coil:	
Receiving coil:	
MR sequence application:	Morpho / Angio / Diffusion / ...
Explored entity:	Anatomical Dataset / fonctionnal / ...
Acquisition contrast:	T1 / T2 / T2* / DP
MR dataset nature:	... weightedMRdataset
MR sequence physics:	... sequence
Dataset name:	...
Dataset comment:	...
MR sequence name:	...
Contrast agent:	NONE / GADO / ...
Contrast agent (manufactured name):	...
Contrast agent concentration:	...
Injected volume:	...
Magnetization transfer:	Yes / No
Parallel acquisition:	Yes / No
Parallel acquisition technique:	GRAPPA / SENSE / ...
Time reduction factor for the in-plane direction:	...
Time reduction factor for the out-of-plane direction:	...
Slice orientation at acquisition:	tra / sag / cor
Axis orientation at acquisition:	bicalleux / ac-pc
Slice order:	center in-out / increase

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Mr Sequence Physics	Sequence name
Hybrid gradient-echo and spin-echo sequence	
Inversion recovery segmented spin-echo sequence	← FLAIR, STIR
Inversion recovery single shot spin-echo sequence	
Inversion recovery single-echo spin-echo sequence	
Magnetization prepared GRE	
Magnetization prepared segmented GRE EPI	
Magnetization prepared single-shot GRE EPI	← Diffusion (GRE)
Magnetization prepared spoiled GRE	← MPRAGE
Multi-echo gradient-echo sequence	← T2 EG
Multi-echo spin-echo sequence	
Refocused GRE	
Segmented GRE EPI	
Segmented spin-echo sequence	
Single shot spin-echo sequence	← Diffusion (SE)
Single-echo gradient-echo sequence	
Single-echo spin-echo sequence	
Single-shot GRE EPI	
Spin echo-echo planar imaging	
Spoiled GRE	
Standard segmented GRE EPI	
Standard segmented spin-echo sequence	← T2 TSE, DP, T1 TSE
Standard single shot spin-echo sequence	
Standard single-echo spin-echo sequence	
Standard single-shot GRE EPI	
Steady state FID-SE	
Steady state GRE FID	
Steady state GRE SE	

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Sequence Name	Acquisition contrast
FLAIR	T2
STIR	T2
MPRAGE	T1
Diffusion (SE)	T2*
Diffusion (GRE)	T2*
T1 TSE	T1
T2 TSE	T2
T2 EG	T2*
DP	DP

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Copyright 2012